

TINAGL1 Antibody, Biotin conjugated



PACO41297

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be implicated in the adrenocortical zonation and in mechanisms for repressing the CYP11B1 gene expression in adrenocortical cells. This is a non catalytic peptidase C1 family protein.

Gene ID:

TINAGL1

Uniprot

Q9GZM7

Synonyms:

Tubulointerstitial nephritis antigen-like (Glucocorticoid-inducible protein 5) (Oxidized LDL-responsive gene 2 protein) (OLRG-2) (Tubulointerstitial nephritis antigen-related protein) (TIN Ag-related protein) (TIN-Ag-RP), TINAGL1, GIS5 LCN7 OLRG2 TINAGL

Immunogen:

Recombinant Human Tubulointerstitial nephritis antigen-like protein (250-467AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A