

# MS4A7 Antibody, Biotin conjugated



PACO41305

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May be involved in signal transduction as a component of a multimeric receptor complex.

**Gene ID:**

MS4A7

**Uniprot**

Q9GZW8

**Synonyms:**

Membrane-spanning 4-domains subfamily A member 7 (CD20 antigen-like 4) (CD20/FC-epsilon-RI-beta family member 4) (Four-span transmembrane protein 2), MS4A7, 4SPAN2 CD20L4 CFFM4

**Immunogen:**

Recombinant Human Membrane-spanning 4-domains subfamily A member 7 protein (1-47AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A