

# UBA5 Antibody, HRP conjugated



PACO41307

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

E1-like enzyme which activates UFM1 and SUMO2.

**Gene ID:**

UBA5

**Uniprot**

Q9GZZ9

**Synonyms:**

Ubiquitin-like modifier-activating enzyme 5 (Ubiquitin-activating enzyme 5) (ThiFP1) (UFM1-activating enzyme) (Ubiquitin-activating enzyme E1 domain-containing protein 1), UBA5, UBE1DC1

**Immunogen:**

Recombinant Human Ubiquitin-like modifier-activating enzyme 5 protein (291-404AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A