

CPVL Antibody, FITC conjugated



PACO41332

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in the digestion of phagocytosed particles in the lysosome, participation in an inflammatory protease cascade, and trimming of peptides for antigen presentation.

Gene ID:

CPVL

Uniprot

Q9H3G5

Synonyms:

Probable serine carboxypeptidase CPVL (EC 3.4.16. -) (Carboxypeptidase, vitellogenic-like) (Vitellogenic carboxypeptidase-like protein) (VCP-like protein) (hVLP), CPVL, VLP

Immunogen:

Recombinant Human Probable serine carboxypeptidase CPVL protein (23-476AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A