

PACO41334

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Can dephosphorylate single and diphosphorylated synthetic MAPK peptides, with preference for the phosphotyrosine and diphosphorylated forms over phosphothreonine.

Gene ID:

DUSP21

Uniprot

Q9H596

Synonyms:

Dual specificity protein phosphatase 21 (EC 3.1.3.16) (EC 3.1.3.48) (Low molecular weight dual specificity phosphatase 21) (LMW-DSP21), DUSP21, LMWDSP21


Immunogen:

Recombinant Human Dual specificity protein phosphatase 21 protein (1-190AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO41334 at dilution of 1:100.