

NRXN3 Antibody, FITC conjugated



PACO41368

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Neuronal cell surface protein that may be involved in cell recognition and cell adhesion. May play a role in angiogenesis.

Gene ID:

NRXN3

Uniprot

Q9HDB5

Synonyms:

Neurexin-3-beta (Neurexin III-beta) [Cleaved into: Neurexin-3-beta, soluble form; Neurexin-3-beta, C-terminal fragment (NRXN3-CTF)], NRXN3, KIAA0743

Immunogen:

Recombinant Human Neurexin-3-β protein (201-432AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A