

ZCCHC17 Antibody, HRP conjugated



PACO41375

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

intermediate filament cytoskeleton, intracellular membrane-bounded organelle, nucleus, poly(A) RNA binding, mitophagy in response to mitochondrial depolarization, positive regulation of defense response to virus by host, xenophagy.

Gene ID:

ZCCHC17

Uniprot

Q9NP64

Synonyms:

Nucleolar protein of 40 kDa (pNO40) (Pnn-interacting nucleolar protein) (Putative S1 RNA-binding domain protein) (PS1D protein) (Zinc finger CCHC domain-containing protein 17), ZCCHC17, PS1D

Immunogen:

Recombinant Human Nucleolar protein of 40 kDa protein (1-241AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A