

# ZCCHC17 Antibody, FITC conjugated



PACO41376

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

intermediate filament cytoskeleton, intracellular membrane-bounded organelle, nucleus, poly(A) RNA binding, mitophagy in response to mitochondrial depolarization, positive regulation of defense response to virus by host, xenophagy.

**Gene ID:**

ZCCHC17

**Uniprot**

Q9NP64

**Synonyms:**

Nucleolar protein of 40 kDa (pNO40) (Pnn-interacting nucleolar protein) (Putative S1 RNA-binding domain protein) (PS1D protein) (Zinc finger CCHC domain-containing protein 17), ZCCHC17, PS1D

**Immunogen:**

Recombinant Human Nucleolar protein of 40 kDa protein (1-241AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A