

EXOSC5 Antibody, HRP conjugated



PACO41391

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

EXOSC5

Uniprot

Q9NQT4

Synonyms:

Exosome complex component RRP46 (Chronic myelogenous leukemia tumor antigen 28) (Exosome component 5) (Ribosomal RNA-processing protein 46) (p12B), EXOSC5, CML28 RRP46

Immunogen:

Recombinant Human Exosome complex component RRP46 protein (1-235AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A