

# EXOSC5 Antibody, Biotin conjugated



PACO41393

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

EXOSC5

**Uniprot**

Q9NQT4

**Synonyms:**

Exosome complex component RRP46 (Chronic myelogenous leukemia tumor antigen 28) (Exosome component 5) (Ribosomal RNA-processing protein 46) (p12B), EXOSC5, CML28 RRP46

**Immunogen:**

Recombinant Human Exosome complex component RRP46 protein (1-235AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A