

# ANKRD10 Antibody



PACO41462

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

regulation of canonical Wnt signaling pathway.

**Gene ID:**

ANKRD10

**Uniprot**

Q9NXR5

**Synonyms:**

Ankyrin repeat domain-containing protein 10, ANKRD10

**Immunogen:**


Recombinant Human Ankyrin repeat domain-containing protein 10 protein (1-220AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Western Blot. Positive WB detected in: Mouse heart tissue, Mouse kidney tissue, Mouse brain tissue. All lanes: ANKRD10 antibody at 4 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 45, 24 kDa. Observed band size: 50 kDa.