

RASL12 Antibody, Biotin conjugated



PACO41489

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RASL12

Uniprot

Q9NYN1

Synonyms:

Ras-like protein family member 12 (Ras-like protein Ris), RASL12, RIS

Immunogen:

Recombinant Human Ras-like protein family member 12 protein (1-266AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A