

# NIPSNAP3B Antibody



PACO41526

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:****Gene ID:**

NIPSNAP3B

**Uniprot**

Q9BS92

**Synonyms:**

Protein NipSnap homolog 3B (NipSnap3B) (SNAP1), NIPSNAP3B, NIPSNAP3

**Immunogen:**


Recombinant Human Protein NipSnap homolog 3B protein (1-247AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO41526 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO41526 at dilution of 1:100.