

DGCR6L Antibody, Biotin conjugated



PACO41557

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in neural crest cell migration into the third and fourth pharyngeal pouches.

Gene ID:

DGCR6L

Uniprot

Q9BY27

Synonyms:

Protein DGCR6L (DiGeorge syndrome critical region 6-like protein), DGCR6L

Immunogen:

Recombinant Human Protein DGCR6L protein (1-220AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A