

SLC25A18 Antibody, FITC conjugated



PACO41600

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in the transport of glutamate across the inner mitochondrial membrane. Glutamate is cotransported with H(+).

Gene ID:

SLC25A18

Uniprot

Q9H1K4

Synonyms:

Mitochondrial glutamate carrier 2 (GC-2) (Glutamate/H(+) symporter 2) (Solute carrier family 25 member 18), SLC25A18, GC2

Immunogen:

Recombinant Human Mitochondrial glutamate carrier 2 protein (246-287AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A