

SH3BGRL3 Antibody, FITC conjugated



PACO41604

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Could act as a modulator of glutaredoxin biological activity.

Gene ID:

SH3BGRL3

Uniprot

Q9H299

Synonyms:

SH3 domain-binding glutamic acid, rich-like protein 3 (SH3 domain-binding protein 1) (SH3BP-1), SH3BGRL3

Immunogen:

Recombinant Human SH3 domain-binding glutamic acid, rich-like protein 3 protein (2-93AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A