

C11orf1 Antibody



PACO41630

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

nucleoplasm, nucleus.

Gene ID:

C11orf1

Uniprot

Q9H5F2

Synonyms:

UPF0686 protein C11orf1, C11orf1


Immunogen:

Recombinant Human UPF0686 protein C11orf1 protein (33-150AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO41630 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO41630 at dilution of 1:100.