

FAM124B Antibody



PACO41642

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

FAM124B

Uniprot

Q9H5Z6

Synonyms:

Protein FAM124B, FAM124B


Immunogen:

Recombinant Human Protein FAM124B protein (146-272AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: FAM124B antibody at 12 μ g/ml. Lane 1: HeLa whole cell lysate. Lane 2: K562 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 51, 32 kDa. Observed band size: 51 kDa.



Immunohistochemistry of paraffin-embedded human lung tissue using PACO41642 at dilution of 1:100.