

# C17orf75 Antibody



PACO41662

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May have a role in spermatogenesis.

**Gene ID:**

C17orf75

**Uniprot**

Q9HAS0

**Synonyms:**

Protein Njmu-R1, C17orf75

**Immunogen:**

Recombinant Human Protein Njmu-R1 protein (201-396AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A