

# TMEM27 Antibody



PACO41666

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200,  
IF:1:50-1:200

**Protein Background:**

Regulator of SNARE complex function. Stimulator of beta cell replication.

**Gene ID:**

TMEM27

**Uniprot**

Q9HBJ8

**Synonyms:**

Collectrin (Transmembrane protein 27), TMEM27

**Immunogen:**


Recombinant Human Collectrin protein (15-141AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images


---



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO41666 at dilution of 1:100.



Immunofluorescent analysis of MCF-7 cells using PACO41666 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO41666 at dilution of 1:100.