

TMEM27 Antibody, HRP conjugated



PACO41667

Product Information

Size:	Protein Background:
50ug	Regulator of SNARE complex function. Stimulator of beta cell replication.
Reactivity:	Gene ID:
Human	TMEM27
Source:	Uniprot
Rabbit	Q9HBJ8
Isotype:	Synonyms:
IgG	Collectrin (Transmembrane protein 27), TMEM27
Applications:	Immunogen:
ELISA	Recombinant Human Collectrin protein (15-141AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A