

TMEM27 Antibody, FITC conjugated



PACO41668

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Regulator of SNARE complex function. Stimulator of beta cell replication.

Gene ID:

TMEM27

Uniprot

Q9HBJ8

Synonyms:

Collectrin (Transmembrane protein 27), TMEM27

Immunogen:

Recombinant Human Collectrin protein (15-141AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A