

VPS45 Antibody, Biotin conjugated



PACO41693

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in vesicle-mediated protein trafficking from the Golgi stack through the trans-Golgi network.

Gene ID:

VPS45

Uniprot

Q9NRW7

Synonyms:

Vacuolar protein sorting-associated protein 45 (h-VPS45) (hVps45), VPS45, VPS45A, VPS45B

Immunogen:

Recombinant Human Vacuolar protein sorting-associated protein 45 protein (101-400AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A