

GPR173 Antibody



PACO41694

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Orphan receptor.

Gene ID:

GPR173

Uniprot

Q9NS66

Synonyms:

Probable G-protein coupled receptor 173 (Super conserved receptor expressed in brain 3), GPR173, SREB3


Immunogen:

Recombinant Human Probable G-protein coupled receptor 173 protein (210-287AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO41694 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO41694 at dilution of 1:100.