

MRPS18A Antibody, HRP conjugated



PACO41707

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

mitochondrial inner membrane, mitochondrial small ribosomal subunit, structural constituent of ribosome, mitochondrial translation, mitochondrial translational elongation, mitochondrial translational initiation, mitochondrial translational termination, organelle organization, translation.

Gene ID:

MRPS18A

Uniprot

Q9NVS2

Synonyms:

39S ribosomal protein S18a, mitochondrial (MRP-S18-a) (Mrps18a) (S18mt-a) (39S ribosomal protein S18-3, mitochondrial) (MRP-S18-3) (Mitochondrial large ribosomal subunit protein bS18a) (Mitochondrial large ribosomal subunit protein mL66), MRPS18A

Immunogen:

Recombinant Human 39S ribosomal protein S18a, mitochondrial protein (1-196AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A