

# TRNAU1AP Antibody, Biotin conjugated



PACO41721

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in the early steps of selenocysteine biosynthesis and tRNA(Sec) charging to the later steps resulting in the cotranslational incorporation of selenocysteine into selenoproteins. Stabilizes the SECISBP2, EEFSEC and tRNA(Sec) complex. May be involved in the methylation of tRNA(Sec). Enhances efficiency of selenoproteins synthesis.

**Gene ID:**

TRNAU1AP

**Uniprot**

Q9NX07

**Synonyms:**

tRNA selenocysteine 1-associated protein 1 (SECp43) (tRNA selenocysteine-associated protein 1), TRNAU1AP, SECP43 TRSPAP1

**Immunogen:**

Recombinant Human tRNA selenocysteine 1-associated protein 1 protein (33-156AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A