

CCDC59 Antibody, FITC conjugated



PACO41764

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of the transcription complexes of the pulmonary surfactant-associated protein-B (SFTPB) and -C (SFTPC). Enhances homeobox protein Nkx-2.1-activated SFTPB and SFTPC promoter activities.

Gene ID:

CCDC59

Uniprot

Q9P031

Synonyms:

Thyroid transcription factor 1-associated protein 26 (TTF-1-associated protein 26) (Coiled-coil domain-containing protein 59) (TTF-1-associated protein BR2), CCDC59, BR22 TAP26

Immunogen:

Recombinant Human Thyroid transcription factor 1-associated protein 26 protein (1-241AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A