

# ASH2L Antibody, HRP conjugated



PACO41771

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Component of the Set1/Ash2 histone methyltransferase (HMT) complex, a complex that specifically methylates Lys-4 of histone H3, but not if the neighboring Lys-9 residue is already methylated. As part of the MLL1/MLL complex it is involved in methylation and dimethylation at Lys-4 of histone H3. May function as a transcriptional regulator. May play a role in hematopoiesis.

**Gene ID:**

ASH2L

**Uniprot**

Q9UBL3

**Synonyms:**

Set1/Ash2 histone methyltransferase complex subunit ASH2 (ASH2-like protein), ASH2L, ASH2L1

**Immunogen:**

Recombinant Human Set1/Ash2 histone methyltransferase complex subunit ASH2 protein (1-534AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A