

PACO41778

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Probable substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex that may function in myogenesis.

Gene ID:

FBXO40

Uniprot

Q9UH90

Synonyms:

F-box only protein 40 (Muscle disease-related protein), FBXO40, FBX40 KIAA1195


Immunogen:

Recombinant Human F-box only protein 40 protein (201-500AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO41778 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human heart tissue using PACO41778 at dilution of 1:100.