

ZDHHC2 Antibody, Biotin conjugated



PACO41809

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Palmitoyltransferase specific for GAP43 and DLG4/PSD95.

Gene ID:

ZDHHC2

Uniprot

Q9UIJ5

Synonyms:

Palmitoyltransferase ZDHHC2 (EC 2.3.1.225) (Reduced expression associated with metastasis protein) (Ream) (Reduced expression in cancer protein) (Rec) (Zinc finger DHHC domain-containing protein 2) (DHHC-2) (Zinc finger protein 372), ZDHHC2, REAM REC ZNF372

Immunogen:

Recombinant Human Palmitoyltransferase ZDHHC2 protein (230-367AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A