

ST6GALNAC2 Antibody, FITC conjugated



PACO41816

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Golgi membrane, sialyltransferase activity, cellular protein metabolic process, dolichol-linked oligosaccharide biosynthetic process, ganglioside biosynthetic process, O-glycan processing, oligosaccharide metabolic process, post-translational protein modification, protein N-linked glycosylation via asparagine, protein O-linked glycosylation.

Gene ID:

ST6GALNAC2

Uniprot

Q9UJ37

Synonyms:

Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 2 (EC 2.4.99. -) (GalNAc alpha-2,6-sialyltransferase II) (ST6GalNAc II) (ST6GalNAcII) (SThM) (Sialyltransferase 7B) (SIAT7-B), ST6GALNAC2, SIAT7B SIATL1 STHM

Immunogen:

Recombinant Human Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 2 protein (101-374AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A