

RABGEF1 Antibody



PACO41818

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Rab effector protein acting as linker between gamma-adaptin, RAB4A or RAB5A. Involved in endocytic membrane fusion and membrane trafficking of recycling endosomes. Stimulates nucleotide exchange on RAB5A. Can act as a ubiquitin ligase.

Gene ID:

RABGEF1

Uniprot

Q9UJ41

Synonyms:

Rab5 GDP/GTP exchange factor (RAP1) (Rabaptin-5-associated exchange factor for Rab5) (Rabex-5), RABGEF1, RABEX5


Immunogen:

Recombinant Human Rab5 GDP/GTP exchange factor protein (1-491AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO41818 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO41818 at dilution of 1:100.