

THTPA Antibody



PACO41830

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Hydrolase highly specific for thiamine triphosphate (ThTP).

Gene ID:

THTPA

Uniprot

Q9BU02

Synonyms:

Thiamine-triphosphatase (ThTPase) (EC 3.6.1.28), THTPA


Immunogen:

Recombinant Human Thiamine-triphosphatase protein (2-230AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO41830 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO41830 at dilution of 1:100.