

CCDC32 Antibody



PACO41838

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Gene ID:

CCDC32

Uniprot

Q9BV29

Synonyms:

Coiled-coil domain-containing protein 32, CCDC32, C15orf57

Immunogen:

Recombinant Human Coiled-coil domain-containing protein 32 protein (1-185AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: CCDC32 antibody at 2 μ g/ml + HL60 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 21, 22, 16 kDa. Observed band size: 21 kDa.