

FNDC11 Antibody



PACO41842

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

extracellular exosome.

Gene ID:

FNDC11

Uniprot

Q9BVV2

Synonyms:

Fibronectin type III domain-containing protein 11, FNDC11, C20orf195


Immunogen:

Recombinant Human Uncharacterized protein C20orf195 protein (10-310AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO41842 at dilution of 1:100.