

TEX13A Antibody



PACO41850

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

TEX13A

Uniprot

Q9BXU3

Synonyms:

Testis-expressed protein 13A, TEX13A


Immunogen:

Recombinant Human Testis-expressed sequence 13A protein (116-380AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO41850 at dilution of 1:100.