

OSBPL6 Antibody



PACO41866

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Weakly binds 25-hydroxycholesterol.

Gene ID:

OSBPL6

Uniprot

Q9BZF3

Synonyms:

Oxysterol-binding protein-related protein 6 (ORP-6) (OSBP-related protein 6), OSBPL6, ORP6

Immunogen:

Recombinant Human Oxysterol-binding protein-related protein 6 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO41866 at dilution of 1:100.