

FAM107B Antibody



PACO41886

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

FAM107B

Uniprot

Q9H098

Synonyms:

Protein FAM107B, FAM107B, C10orf45


Immunogen:

Recombinant Human Protein FAM107B protein (2-131AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: FAM107B antibody at 8 μ g/ml + Mouse heart tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 16, 35 kDa. Observed band size: 16 kDa.