

# CRISPLD2 Antibody, FITC conjugated



PACO41892

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Promotes matrix assembly.

**Gene ID:**

CRISPLD2

**Uniprot**

Q9H0B8

**Synonyms:**

Cysteine-rich secretory protein LCCL domain-containing 2 (Cysteine-rich secretory protein 11) (CRISP-11) (LCCL domain-containing cysteine-rich secretory protein 2), CRISPLD2, CRISP11 LCRISP2

**Immunogen:**

Recombinant Human Cysteine-rich secretory protein LCCL domain-containing 2 protein (40-305AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A