

# OSBPL2 Antibody, HRP conjugated



PACO41903

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Binds phospholipids; exhibits strong binding to phosphatidic acid, and weak binding to phosphatidylinositol 3-phosphate. Binds 25-hydroxycholesterol.

**Gene ID:**

OSBPL2

**Uniprot**

Q9H1P3

**Synonyms:**

Oxysterol-binding protein-related protein 2 (ORP-2) (OSBP-related protein 2), OSBPL2, KIAA0772 ORP2

**Immunogen:**

Recombinant Human Oxysterol-binding protein-related protein 2 protein (251-480AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A