

CRISPLD1 Antibody, Biotin conjugated



PACO41913

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

extracellular exosome, face morphogenesis.

Gene ID:

CRISPLD1

Uniprot

Q9H336

Synonyms:

Cysteine-rich secretory protein LCCL domain-containing 1 (CocoaCrisp) (Cysteine-rich secretory protein 10) (CRISP-10) (LCCL domain-containing cysteine-rich secretory protein 1) (Trypsin inhibitor HI), CRISPLD1, CRISP10 LCRISP1

Immunogen:

Recombinant Human Cysteine-rich secretory protein LCCL domain-containing 1 protein (24-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A