

CUEDC2 Antibody, Biotin conjugated



PACO41917

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Down-regulates ESR1 protein levels through the ubiquitination-proteasome pathway, regardless of the presence of 17 beta-estradiol. Also involved in 17 beta-estradiol-induced ESR1 degradation. Controls PGR protein levels through a similar mechanism.

Gene ID:

CUEDC2

Uniprot

Q9H467

Synonyms:

CUE domain-containing protein 2, CUEDC2, C10orf66

Immunogen:

Recombinant Human CUE domain-containing protein 2 protein (1-226AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A