

# PCIF1 Antibody, HRP conjugated



PACO41923

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play a role in transcription elongation or in coupling transcription to pre-mRNA processing through its association with the phosphorylated C-terminal domain (CTD) of RNAPII largest subunit.

**Gene ID:**

PCIF1

**Uniprot**

Q9H4Z3

**Synonyms:**

Phosphorylated CTD-interacting factor 1, PCIF1, C20orf67

**Immunogen:**

Recombinant Human Phosphorylated CTD-interacting factor 1 protein (201-500AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A