

# NMUR1 Antibody, HRP conjugated



PACO41943

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Receptor for the neuromedin-U and neuromedin-S neuropeptides.

**Gene ID:**

NMUR1

**Uniprot**

Q9HB89

**Synonyms:**

Neuromedin-U receptor 1 (NMU-R1) (G-protein coupled receptor 66) (G-protein coupled receptor FM-3), NMUR1, GPR66

**Immunogen:**

Recombinant Human Neuromedin-U receptor 1 protein (360-426AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A