

Flii Antibody, HRP conjugated



PACO41959

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a role as coactivator in transcriptional activation by hormone-activated nuclear receptors (NR) and acts in cooperation with NCOA2 and CARM1. Involved in estrogen hormone signaling. By similarity. Essential for early embryonic development. May play a role in regulation of cytoskeletal rearrangements involved in cytokinesis and cell migration, by inhibiting Rac1-dependent paxillin phosphorylation.

Gene ID:

Flii

Uniprot

Q9JJ28

Synonyms:

Protein flightless-1 homolog, Flii, Fli1 Fliih

Immunogen:

Recombinant Mouse Protein flightless-1 homolog protein (495-827AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A