

# ZNF71 Antibody



PACO41962

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZNF71

**Uniprot**

Q9NQZ8

**Synonyms:**

Endothelial zinc finger protein induced by tumor necrosis factor alpha (Zinc finger protein 71), ZNF71, EZFIT

**Immunogen:**


Recombinant Human Endothelial zinc finger protein induced by tumor necrosis factor &alpha; protein (101-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO41962 at dilution of 1:100.