

PNO1 Antibody, Biotin conjugated



PACO41969

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

nucleolus, nucleus, poly(A) RNA binding.

Gene ID:

PNO1

Uniprot

Q9NRX1

Synonyms:

RNA-binding protein PNO1, PNO1

Immunogen:

Recombinant Human RNA-binding protein PNO1 protein (1-252AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A