

# ASNSD1 Antibody



PACO41990

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

cytosol, asparagine synthase (glutamine-hydrolyzing) activity, protein homodimerization activity, asparagine biosynthetic process, glutamine metabolic process.

**Gene ID:**

ASNSD1

**Uniprot**

Q9NWL6

**Synonyms:**

Asparagine synthetase domain-containing protein 1 (HCV NS3-transactivated protein 1), ASNSD1, NS3TP1

**Immunogen:**


Recombinant Human Asparagine synthetase domain-containing protein 1 protein (258-357AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human liver cancer using PACO41990 at dilution of 1:100.