

ASNSD1 Antibody, FITC conjugated



PACO41992

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytosol, asparagine synthase (glutamine-hydrolyzing) activity, protein homodimerization activity, asparagine biosynthetic process, glutamine metabolic process.

Gene ID:

ASNSD1

Uniprot

Q9NWL6

Synonyms:

Asparagine synthetase domain-containing protein 1 (HCV NS3-transactivated protein 1), ASNSD1, NS3TP1

Immunogen:

Recombinant Human Asparagine synthetase domain-containing protein 1 protein (258-357AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A