

MRPL36 Antibody



PACO42010

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of the large subunit of the mitochondrial ribosome.

Gene ID:

MRPL36

Uniprot

Q9P0J6

Synonyms:

39S ribosomal protein L36, mitochondrial (L36mt) (MRP-L36) (BRCA1-interacting protein 1) (Mitochondrial large ribosomal subunit protein bL36m), MRPL36, BRIP1

Immunogen:

Recombinant Human 39S ribosomal protein L36, mitochondrial protein (2-103AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A